

Infinium™ Omni2.5Exome-8 v1.4 BeadChip

Comprehensive coverage and functional exonic content for genotyping, GWAS, and CNV analysis.

Overview

The Infinium Omni2.5Exome-8 v1.4 BeadChip (Figure 1) delivers comprehensive coverage of common, rare, and exonic single nucleotide polymorphism (SNP) content from the 1000 Genomes Project,¹ providing maximum genomic information of diverse world populations (Table 1 and Table 2). With combined markers from the Infinium Omni2.5-8 and Infinium Exome-24 BeadChips, the Infinium Omni2.5Exome-8 v1.4 BeadChip is a powerful tool for next-generation genotyping and genome-wide association studies (GWAS) (Table 3). Using the proven HiScan™ or iScan™ System and integrated analysis software, the 8-sample Infinium Omni2.5Exome-8 v1.4 BeadChip offers high throughput, optimized tag SNPs, functional exonic content, and fully supported copy-number variation (CNV) analysis. Combined with convenient packaging and a streamlined PCR-free protocol, the Infinium Omni2.5Exome-8 v1.4 BeadChip provides a comprehensive DNA analysis solution.



Figure 1: The Infinium Omni2.5Exome-8 v1.4 BeadChip— The Infinium Omni2.5Exome-8 v1.4 BeadChip supports rapid, cost-effective GWAS studies with coverage of common, rare, and functional variants from the 1000 Genomes Project.

Table 1: Product information

Feature	Description		
Species	Human		
Total number of markers	2,619,927		
Number of samples per BeadChip	8 Samples		
DNA input requirement	200 ng		
Assay chemistry	Infinium LCG		
Instrument support	iScan or HiScan System		
Sample throughput ^a	~ 1067 samples/week		
Scan time per sample	iScan System	HiScan System	
	11.3 min	6.5 min	
Data performance	Value ^b	Product Specification ^d	
Call rate	99.8%	> 99% avg.	
Reproducibility	99.99%	> 99.9%	
Log R deviation	0.11	< 0.30 ^c	
Spacing			
Spacing (kb)	Mean	Median	90th% ^c
	1.11	0.56	2.63

- Estimate assumes 1 iScan System, 1 AutoLoader 2.x, 2 Tecan robots, and a 5-day work week.
- Values are derived from genotyping 458 HapMap reference samples.
- Value expected for typical projects using standard Illumina protocols. Tumor samples and samples prepared by methods other than standard Illumina protocols are excluded.
- Excludes Y chromosome markers for female samples.

Table 2: LD $r^2 \geq 0.80$ from 1000G^a at various MAF thresholds

1000G population ^b	LD coverage ($r^2 \geq 0.80$)		
	MAF $\geq 1\%$	MAF $\geq 2.5\%$	MAF $\geq 5\%$
AFR	0.54	0.64	0.69
AMR	0.74	0.82	0.86
EAS	0.79	0.86	0.89
EUR	0.79	0.86	0.89
SAS	0.75	0.83	0.87

- Compared against Phase 3, version 5 of the 1000 Genomes Project (1000G). www.1000genomes.org. Accessed July 2016.
- See www.1000genomes.org/category/frequently-asked-questions/population

Abbreviations: LD, linkage disequilibrium; MAF, minor allele frequency; AFR, African; AMR, Ad Mixed American; EAS, East Asian; EUR, European; SAS, South Asian.

Table 3: Marker information

Marker categories	No. of markers		
RefSeq ^a genes	1,337,616		
RefSeq exons	363,435		
RefSeq promoter regions	57,433		
ADME ^b genes	23,745		
Extended MHC	22,480		
COSMIC ^c genes	1,190,730		
Genes in Gene Ontology ^d	309,642		
Nonsense markers ^e	5506		
Missense markers ^e	237,687		
Synonymous markers ^e	41,606		
Silent markers ^f	78,673		
Mitochondrial markers ^f	400		
Indels ^f	172		
Sex chromosomes ^f	X	Y	PAR/homologous
	58,301	2464	2710

- a. RefSeq - NCBI Reference Sequence Database. www.ncbi.nlm.nih.gov/refseq. Accessed September 2016.
- b. PharmaADME Gene List. www.pharmaadme.org. Accessed August 2014.
- c. Catalog of somatic mutations in cancer. cancer.sanger.uk/cosmic. Accessed July 2016.
- d. Gene Ontology Consortium. www.geneontology.org. Accessed July 2016.
- e. Compared against the University of California, Santa Cruz (UCSC) Genome Browser. genome.ucsc.edu. Accessed August 2014.
- f. NCBI Genome Reference Consortium, Version GRCh37. www.ncbi.nlm.nih.gov/grc/human. Accessed July 2016.

Abbreviations: indel: insertion/deletion; PAR: pseudoautosomal region.

Ordering information

Infinium Omni2.5Exome-8 v1.4 Kit	Catalog no.
16 samples	20024552
48 samples	20024553
96 samples	20024554
384 samples	20024555

Learn more

To learn more about the Infinium Omni2.5Exome-8 v1.4 BeadChip and other Illumina genotyping products and services, visit www.illumina.com/genotyping

References

- 1000 Genomes Project, www.1000genomes.org. Accessed April 2014.