

Amplify DNA

- 1 Add DNA into either of the following to create a DNA plate:
 - ▶ Midi plate: 20 µl to each DNA well
 - ▶ TCY plate: 10 µl to each DNA well
- 2 Select **MSA6 Tasks | Make MSA6**.
- 3 Select the DNA plate type.
- 4 Enter the **Number of DNA plates**.
- 5 Place the MA1, RPM, and MSM tubes in the robot tube rack.
- 6 Pour 15 ml NaOH into a trough and place on the robot bed.
- 7 Place DNA and MSA6 plates on robot bed.
- 8 Select **Run**.
- 9 Enter the barcode of each DNA plate.
- 10 Place the DNA plates on the robot bed and select **OK**.
- 11 Vortex the MSA6 plate at 1600 rpm for 1 minute.
- 12 Centrifuge at 280 × g.
- 13 Remove the cap mat, place the MSA6 plate on the robot bed, and select **OK**.
- 14 When complete, select **OK**.
- 15 Remove and seal the MSA6 plate.
- 16 Centrifuge at 280 × g.

Incubate DNA

- 1 **[LIMS]** Select **Infinium LCG**
 - a Scan the barcodes.
- 2 Incubate the MSA6 plate for 20–24 hours at 37°C.

Fragment DNA

- 1 Pulse centrifuge the MSA6 plate at 280 × g.
- 2 Select **MSA6 Tasks | Fragment MSA6**.
- 3 Place the MSA6 plate on the robot bed.
- 4 Place FMS tubes in the robot tube rack.
- 5 Select **Run**.
- 6 When complete, select **OK**.
- 7 Remove the plate and seal with a cap mat.
- 8 Vortex at 1600 rpm for 1 minute.
- 9 Pulse centrifuge at 280 × g.
- 10 Place on the 37°C heat block for 1 hour.

SAFE STOPPING POINT

If you are stopping, seal the plate(s), and store at -25°C to -15°C.

Precipitate DNA

- 1 Select **MSA6 Tasks | Precip MSA6**.
- 2 Pulse centrifuge the sealed plate at 280 × g.
- 3 Place the MSA6 plate on the robot bed.
- 4 Place a half reservoir in the frame, and add PM1 as follows:
 - ▶ For 48 samples, add 1 tube PM1
 - ▶ For 96 samples, add 2 tubes PM1
- 5 Place a full reservoir in the frame, and add 2-propanol as follows:
 - ▶ For 48 samples, add 20 ml 2-propanol
 - ▶ For 96 samples, add 40 ml 2-propanol
- 6 Select **Run**.
- 7 Remove the MSA6 plate from the robot bed. Do not select **OK**.
- 8 Vortex at 1600 rpm for 1 minute.
- 9 Incubate on the heat block for 5 minutes.
- 10 Centrifuge at 280 × g for 1 minute.
- 11 Set the centrifuge at 4°C.
- 12 Place the MSA6 plate on the robot bed.
- 13 Select **OK**.
- 14 Remove the MSA6 plate from the robot bed and seal.
- 15 Invert 10 times to mix.
- 16 Incubate at 4°C for 30 minutes.
- 17 Place in the centrifuge.
- 18 Centrifuge at 3000 × g for 20 minutes.
- 19 Remove MSA6 plate.
- 20 Make sure that a blue pellet is present.
- 21 Remove and discard the cap mat.
- 22 Quickly invert the plate and drain the supernatant.
- 23 Firmly tap until all wells are free of liquid.

- 24 Place the plate on a tube rack for 1 hour at room temperature.
- 25 Make sure that a blue pellet is still present.

SAFE STOPPING POINT

If you are stopping, seal the plate(s), and store at -25°C to -15°C.

Resuspend DNA

- 1 Select **MSA6 Tasks | Resuspend MSA6**.
- 2 Place the MSA6 plate on the robot bed.
- 3 Place a quarter reservoir in the frame, and add RA1 as follows:
 - ▶ For 48 samples, add 4.5 ml RA1
 - ▶ For 96 samples, add 9 ml RA1
- 4 Select **Run**.
- 5 Remove the MSA6 plate from the robot deck.
- 6 Apply a foil seal to the MSA6 plate.
- 7 Incubate in the Illumina Hybridization Oven for 1 hour.
- 8 Vortex at 1800 rpm for 1 minute.
- 9 Make sure that the pellets are resuspended.
- 10 Pulse centrifuge at 280 × g.

SAFE STOPPING POINT

If you are stopping, store sealed MSA6 plate(s) at 2°C to 8°C for up to 24 hours. If more than 24 hours, store at -25°C to -15°C.

Store sealed RA1 at -25°C to -15°C. If RA1 will be used the next day, seal it, and store it overnight at 4°C.

Hybridize DNA to the BeadChip

- 1 Incubate the MSA6 plate on the heat block for 20 minutes.
- 2 Cool at room temperature for 30 minutes.
- 3 Pulse centrifuge at 280 × g.
- 4 Place the gasket into the hybridization chamber.
- 5 Add 400 µl PB2 into each reservoir.
- 6 Place the hybridization chamber insert into the hybridization chamber.
- 7 Immediately cover the chamber with the lid.
- 8 **[LIMS]** Select **Select Infinium LCG | Confirm for Hyb.**
- 9 **[LIMS]** Scan the barcodes.
- 10 Remove all BeadChips from packaging.
- 11 Place BeadChips into the robot BeadChip alignment fixtures.
- 12 Select **MSA6 Tasks Tasks | Hyb.**
 - a Select the 24-sample BeadChip.
 - b Enter the **Number of MSA6 plates.**
- 13 Place the robot BeadChip alignment fixtures onto the robot deck.
- 14 Pulse centrifuge the MSA6 plate at 280 × g.
- 15 Place the MSA6 plate onto the robot deck.
- 16 Select **Run.**
- 17 Place each robot tip alignment guide on top of each robot BeadChip alignment fixture.
- 18 To start the run, select **OK.**
- 19 When complete, select **OK.**
- 20 Remove the robot BeadChip alignment fixtures.
- 21 Place each BeadChip in a hybridization chamber insert.
- 22 Place the lid on the chamber and secure with the metal clamps.

- 23 **[LIMS]** Select **Infinium LCG | Prepare Hyb Chamber.**
 - a Scan the barcodes.
- 24 Incubate at 48°C for 16–24 hours.

Prepare for Next Day

- 1 Add 330 ml fresh 100% EtOH to the XC4 bottle.
- 2 Vigorously shake to resuspend.
- 3 Leave the bottle upright on the lab bench overnight.
- 4 Soak the robot tip alignment guides in 1% aqueous Alconox solution.
- 5 Rinse and dry the robot tip alignment guides.

Wash BeadChips

- 1 Submerge the wash rack in the PB1 wash.
- 2 Remove the hybridization insert.
- 3 Remove the BeadChips.
- 4 Remove the cover seals from the BeadChips.
- 5 Place the BeadChips into the submerged wash rack.
- 6 Move the wash rack up and down for 1 minute.
- 7 Move the wash rack to the next PB1 Wash.
- 8 Move the wash rack up and down for 1 minute.
- 9 Confirm that you are using the correct Infinium LCG glass back plates and spacers.
- 10 Fill the BeadChip alignment fixture with 150 ml PB1.
- 11 For each BeadChip, place one black frame into the BeadChip alignment fixture.
- 12 Place each BeadChip into a black frame.
- 13 Place a **clear** spacer onto the top of each BeadChip.
- 14 Place the alignment bar onto the alignment fixture.
- 15 Place a clean glass back plate on top of each clear spacer.
- 16 Secure each flow-through chamber assembly with metal clamps.
- 17 Remove the assembled flow-through chamber from the alignment fixture.
- 18 Trim the spacers from each end of the assembly.
- 19 Leave assembled flow-through chambers on the lab bench.

- 20 Wash the hybridization chamber reservoirs with DI H₂O.

Extend and Stain BeadChips

- 1 Fill the water circulator.
- 2 Select **Robot QC Tasks | Circulator Manager** to set to 44°C.
- 3 Select **XStain Tasks | XStain LCG BeadChip**.
- 4 Add the following reagents to reservoirs:

Reagent	# BeadChips	Volume
95% formamide/1 mM EDTA	1–8	15 ml
	9–16	17 ml
	17–24	25 ml
RA1	1–8	10 ml
	9–16	20 ml
	17–24	30 ml
XC3	1–8	50 ml
	9–16	100 ml
	17–24	150 ml

- 5 Invert the LX1, LX2, EML, SML, and ATM tubes to mix. Remove the caps, and place on the robot deck.
- 6 Enter the number of BeadChips.
- 7 Select **Run**.
- 8 [Non-LIMS]Enter the stain temperature listed on the SML tube.
- 9 Place the flow-through chambers into the chamber rack.
- 10 Select **OK**.
- 11 Remove the flow-through chambers from the chamber rack.
- 12 Set up two top-loading wash dishes labeled PB1 and XC4.

- 13 Add 310 ml PB1 to the PB1 wash dish.
- 14 Submerge the staining rack in the wash dish.
- 15 Leave the staining rack in the wash dish.
- 16 Disassemble each flow-through chamber.
- 17 Place the BeadChips into the submerged staining rack.
- 18 Slowly lift the staining rack 10 times.
- 19 Soak for 5 minutes.
- 20 Vigorously shake the XC4 bottle.
- 21 Add 310 ml XC4 to the XC4 wash dish and cover.
- 22 Transfer the staining rack from the PB1 to the XC4.
- 23 Slowly lift the staining rack 10 times.
- 24 Soak for 5 minutes.
- 25 Remove the staining rack and place it onto the tube rack.
- 26 Dry each BeadChip as follows.
 - a Grip the BeadChip by the barcode end.
 - b Place onto a tube rack with the barcode facing up and toward you.
- 27 Place the tube rack into the vacuum desiccator.
- 28 Dry the BeadChips for 50–55 minutes at 675 mm Hg (0.9 bar).
- 29 **[LIMS]** Select **Infinium LCG | Coat**.
 - a Scan the barcodes.