

## Tagment Genomic DNA

- 1 Quantify gDNA using a fluorometric method.
- 2 Dilute gDNA in Tris-HCl 10 mM, pH 8.5 to a final volume of 10  $\mu$ l at 5 ng/ $\mu$ l.
- 3 Add the following to a new plate.
  - ▶ Normalized gDNA (10  $\mu$ l)
  - ▶ TD (25  $\mu$ l)
  - ▶ TDE1 (15  $\mu$ l)
- 4 Shake at 1800 rpm for 1 minute.
- 5 Centrifuge at 280  $\times$  g for 1 minute.
- 6 Place on the 58°C microheating system with the lid closed for 10 minutes.
- 7 Add 15  $\mu$ l ST.
- 8 Shake at 1800 rpm for 1 minute.
- 9 Centrifuge at 280  $\times$  g for 1 minute.
- 10 Incubate at room temperature for 4 minutes.

## Clean Up Tagmented DNA

- 1 Add 65  $\mu$ l SPB.
- 2 Shake at 1800 rpm for 1 minute.
- 3 Incubate at room temperature for 8 minutes.
- 4 Centrifuge at 280  $\times$  g for 1 minute.
- 5 Place on a magnetic stand until liquid is clear.
- 6 Remove and discard all supernatant.
- 7 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 8 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 9 Air-dry on the magnetic stand for 10 minutes.
- 10 Remove from the magnetic stand.
- 11 Add 22.5  $\mu$ l RSB.
- 12 Shake at 1800 rpm for 1 minute.
- 13 Incubate at room temperature for 2 minutes.
- 14 Centrifuge at 280  $\times$  g for 1 minute.
- 15 Place on a magnetic stand until liquid is clear.
- 16 Transfer 20  $\mu$ l supernatant.

## Amplify Tagmented DNA

- 1 Arrange Index 1 (i7) adapters in columns 1–12.
- 2 Arrange Index 2 (i5) adapters in rows A–H.
- 3 Place the plate on the TruSeq Index Plate Fixture.
- 4 Add 5  $\mu$ l of each Index 1 adapter down each column.
- 5 Add 5  $\mu$ l of each Index 2 adapter across each row.
- 6 Add 20  $\mu$ l NLM.
- 7 Shake at 1200 rpm for 1 minute.
- 8 Centrifuge at 280  $\times$  g for 1 minute.
- 9 Place on the thermal cycler and run the NLM AMP program.

**SAFE STOPPING POINT**

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

## Clean Up Amplified DNA

- 1 Centrifuge at 280 × g for 1 minute.
- 2 Transfer 50 µl supernatant.
- 3 Add 90 µl SPB.
- 4 Shake at 1800 rpm for 1 minute.
- 5 Incubate at room temperature for 10 minutes.
- 6 Centrifuge at 280 × g for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Wash 2 times with 200 µl 80% EtOH.
- 10 Use a 20 µl pipette to remove residual EtOH.
- 11 Air-dry on the magnetic stand for 10 minutes.
- 12 Add 27 µl RSB.
- 13 Shake at 1800 rpm for 1 minute.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at 280 × g for 1 minute.
- 16 Place on a magnetic stand until liquid is clear.
- 17 Transfer 25 µl supernatant.
- 18 Quantify the library using a fluorometric method.

## SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 14 days.

## Hybridize Probes

- 1 Combine 500 ng of each DNA library. Make sure that each library has a unique index.
  - ▶ For total volume > 40 µl, concentrate the pooled sample to 40 µl.
  - ▶ For total volume < 40 µl, increase the volume to 40 µl with RSB.
- 2 Add the following to a new plate.
  - ▶ DNA library sample or pool (40 µl)
  - ▶ EHB (50 µl)
  - ▶ TCO (10 µl)
- 3 Shake at 1200 rpm for 1 minute.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place on the thermal cycler and run the NRC HYB program.
- 6 Keep at the 58°C holding temperature for at least 90 minutes and up to 24 hours.

## Capture Hybridized Probes

- 1 Centrifuge at 280 × g for 1 minute.
- 2 Transfer all volumes.
- 3 Add 250 µl SMB.
- 4 Shake at 1200 rpm for 5 minutes.
- 5 Incubate at room temperature for 25 minutes.
- 6 Centrifuge at 280 × g for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Remove from the magnetic stand.
- 10 Wash 2 times with 200 µl EWS.
- 11 Mix 28.5 µl EE1 and 1.5 µl HP3, and then vortex.
- 12 Add 23 µl elution premix.
- 13 Shake at 1800 rpm for 2 minutes.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at 280 × g for 1 minute.
- 16 Place on a magnetic stand until liquid is clear.
- 17 Transfer 21 µl supernatant.
- 18 Add 4 µl ET2.
- 19 Shake at 1200 rpm for 1 minute.
- 20 Centrifuge at 280 × g for 1 minute.

## SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

## Perform Second Hybridization

- 1 Add the following.
  - ▶ RSB (15  $\mu$ l)
  - ▶ EHB (50  $\mu$ l)
  - ▶ TCO (10  $\mu$ l)
- 2 Shake at 1200 rpm for 1 minute.
- 3 Centrifuge at 280  $\times$  g for 1 minute.
- 4 Place on the thermal cycler and run the NRC HYB program.
- 5 Keep at the 58°C holding temperature for at least 14.5 hours and up to 24 hours.

## Perform Second Capture

- 1 Centrifuge at 280  $\times$  g for 1 minute.
- 2 Transfer supernatant.
- 3 Add 250  $\mu$ l SMB.
- 4 Shake at 1200 rpm for 5 minutes.
- 5 Incubate at room temperature for 25 minutes.
- 6 Centrifuge at 280  $\times$  g for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Remove from the magnetic stand.
- 10 Wash 2 times with 200  $\mu$ l EWS.
- 11 Mix 28.5  $\mu$ l EE1 and 1.5  $\mu$ l HP3, and then vortex.
- 12 Add 23  $\mu$ l elution premix.
- 13 Shake at 1800 rpm for 2 minutes.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at 280  $\times$  g for 1 minute.
- 16 Place on a magnetic stand until liquid is clear.
- 17 Transfer 21  $\mu$ l supernatant.
- 18 Add 4  $\mu$ l ET2.
- 19 Shake at 1800 rpm for 1 minute.
- 20 Centrifuge at 280  $\times$  g for 1 minute.

## Clean Up Captured Library

- 1 Add 45  $\mu$ l SPB.
- 2 Shake at 1800 rpm for 1 minute.
- 3 Incubate at room temperature for 10 minutes.
- 4 Centrifuge at 280  $\times$  g for 1 minute.
- 5 Place on a magnetic stand until liquid is clear.
- 6 Remove and discard all supernatant.
- 7 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 8 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 9 Air-dry for 10 minutes.
- 10 Add 27.5  $\mu$ l RSB.
- 11 Shake at 1800 rpm for 1 minute.
- 12 Incubate at room temperature for 2 minutes.
- 13 Centrifuge at 280  $\times$  g for 1 minute.
- 14 Place on a magnetic stand until liquid is clear.
- 15 Transfer 25  $\mu$ l supernatant.

**SAFE STOPPING POINT**

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

### Amplify Enriched Library

- 1 Add 5 µl PPC.
- 2 Add 20 µl NEM.
- 3 Shake at 1200 rpm for 1 minute.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place on the thermal cycler and run the NEM AMP12 program.

#### SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days.

### Clean Up Amplified Enriched Library

- 1 Centrifuge at 280 × g for 1 minute.
- 2 Transfer 50 µl.
- 3 Add 90 µl SPB.
- 4 Shake at 1800 rpm for 1 minute.
- 5 Incubate at room temperature for 10 minutes.
- 6 Centrifuge at 280 × g for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Wash 2 times with 200 µl 80% EtOH.
- 10 Use a 20 µl pipette to remove residual EtOH.
- 11 Air-dry on the magnetic stand for 10 minutes.
- 12 Add 32 µl RSB.
- 13 Shake at 1800 rpm for 1 minute.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at 280 × g for 1 minute.
- 16 Place on a magnetic stand until liquid is clear.
- 17 Transfer 30 µl supernatant.

#### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

### Check Enriched Libraries

- 1 Quantify using a fluorometric method.
- 2 If the concentration is higher than the quantitative range for the High Sensitivity DNA chip, dilute the library 1:10 with RSB.
- 3 Run 1 µl diluted using a High Sensitivity DNA chip.
- 4 Denature and dilute pooled libraries to the loading concentration for the instrument you are using. See the denature and dilute libraries guide for your instrument.

## Acronyms

| Acronym | Definition                          |
|---------|-------------------------------------|
| EE1     | Enrichment Elution Buffer 1         |
| EHB     | Enrichment Hybridization Buffer     |
| ET2     | Elute Target Buffer 2               |
| EWS     | Enrichment Wash Solution            |
| NEC1    | Nextera Enriched Clean Up Plate 1   |
| NEC2    | Nextera Enriched Clean Up Plate 2   |
| NEH1    | Nextera Enrichment Hyb Plate 1      |
| NEH2    | Nextera Enrichment Hyb Plate 2      |
| NEL     | Nextera Enrichment Library Plate    |
| NEM     | Enrichment Amp Mix                  |
| NEW1    | Nextera Enrichment Wash Plate 1     |
| NEW2    | Nextera Enrichment Wash Plate 2     |
| NIL     | Nextera Index Library Plate         |
| NLA     | Nextera Library Amplification Plate |
| NLC     | Nextera Library Clean Up Plate      |
| NLM     | Library Amp Mix                     |
| NLT     | Nextera Library Tagment Plate       |
| PPC     | PCR Primer Cocktail                 |
| RSB     | Resuspension Buffer                 |
| SMB     | Streptavidin Magnetic Beads         |

| Acronym | Definition                |
|---------|---------------------------|
| SPB     | Sample Purification Beads |
| ST      | Stop Tagment Buffer       |
| TCO     | TruSight Cardio Oligos    |
| TD      | Tagment DNA Buffer        |
| TDE1    | Tagment DNA Enzyme TDE    |