

Local Run Manager 16S Metagenomics Analysis Module v2

Workflow Guide

This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY, AND WILL VOID ANY WARRANTY APPLICABLE TO THE PRODUCT(S).

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2022 Illumina, Inc. All rights reserved.

All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

Table of Contents

Overview	1
Compatible Library Types	1
About This Guide	1
Set Parameters	2
Set Run Parameters	2
Specify Run Settings	2
Specify Module-Specific Settings	2
Specify Samples for the Run	3
Analysis Methods	5
Demultiplexing	5
FASTQ File Generation	5
Classification of Reads	5
Statistics Reporting	6
View Analysis Results	7
Analysis Report	8
Sample Configuration	8
Sample Information	8
Classification Results	8
Analysis Output Files	10
Supplementary Output Files	10
Custom Analysis Settings	11
Add a Custom Analysis Setting	11
Sample Sheet Settings for Analysis	11
Technical Assistance	12

Overview

The Local Run Manager 16S Metagenomics analysis module classifies organisms from a metagenomic sample by amplifying specific regions in the 16S ribosomal RNA. This workflow is exclusive to Prokaryotes, which includes Bacteria and Archaea. The 16S Metagenomics workflow generates a classification of reads at several taxonomic levels: kingdom, phylum, class, order, family, and genus or species.

The 16S Metagenomics workflow demultiplexes indexed reads, generates FASTQ files, and then classifies reads.

The Local Run Manager 16S Metagenomics (or later) analysis module can only be run on Local Run Manager v3.0 (or later). The analysis module requires the library prep and index kits that are used in the run to be decoupled.

Compatible Library Types

The 16S Metagenomics analysis module is compatible with specific library types represented by library kit categories on the Create Run screen. For a current list of compatible library kits, see the [Platform software Name support page](#) on the Illumina website.

About This Guide

This guide provides instructions for setting up run parameters for sequencing and analysis parameters for the 16S Metagenomics analysis module. For information about the Platform software Name dashboard and system settings, see the Local Run Manager v3 Software Guide (document #100000011492).

Set Parameters

Set Run Parameters

1. If needed, log in to Local Run Manager.
2. Select **Create Run**, and select **16S Metagenomics**.
3. Enter a run name that identifies the run from sequencing through analysis.
The run name must be unique, and cannot be a name used for a previous run set up. The run name can contain up to 40 alphanumeric characters, spaces, and the following special characters:
`~!@#\$\$%-_{}`.
4. [Optional] Enter a run description to identify the run.
The run description can contain up to 150 alphanumeric characters, spaces, and the following special characters: `~!@#\$\$%-_{}`.

Specify Run Settings

1. Select the library prep kit from the Library Prep Kit drop-down list.
2. Select the index kit from the Index Kit drop-down list.
3. Specify the number of index reads.
 - **0** for a run with no indexing
 - **1** for a single-indexed run
 - **2** for a dual-indexed runUnsupported index reads for your index kit are automatically disabled.
4. Specify a read type: **Single Read** or **Paired End**.
If your index kit supports only one option, the read type is automatically selected.
5. Enter the number of cycles for the run.
6. [Optional] If using custom primers, specify their information.
Custom primer options may vary based on your instrument or Local Run Manager implementation.

Specify Module-Specific Settings

1. If applicable, set the Reverse Complement option.
Enabling this setting makes all reads reverse-complemented as they are written to FASTQ files. By default, the reverse complement setting is turned off.

2. Select the taxonomy database to use for the run.
 - **Species-level Classification**—Use an Illumina-developed species-level classification taxonomy database.
 - **Custom Taxonomy Database**—Use a custom taxonomy database for taxonomic classification.
3. If you selected **Custom Taxonomy Database** as the taxonomy database, enter the full file path for the custom taxonomy database file.

Specify Samples for the Run

Specify samples for the run using the following options:

- **Enter samples manually**—Use the blank table at the bottom of the Create Run screen.
- **Import sample sheet**—Navigate to an external file in a comma-separated values (*.csv) format.

After you have populated the samples table, you can export the sample information to an external file. You can use this file as a reference when preparing libraries or import the file when configuring another run.

Enter Samples Manually

To enter sample information manually, you must first select a Library Prep Kit and index kit in the Run Settings section.

1. Adjust the samples table to an appropriate number of rows.
 - In the Rows field, use the up/down arrows or enter a number to specify the number of rows to add to the table. Select the Add Row icon to add the rows to the table.
 - Select the Delete Row icon to delete a row.
 - Right-click on a row in the table and use the commands in the contextual menu.
2. Enter a unique sample ID in the Sample ID field.

Only use alphanumeric characters or dashes. The Sample ID field has a 40 character limit. Do not include spaces. Do not include an underscore as it will be converted to dashes in fastq file names.
3. Enter a sample name.

This field is required when a sample sheet with a Sample_Name column is imported.
4. [Optional] Enter a sample description in the Sample Description field.

Only use alphanumeric characters, dashes, or underscores. Do not include spaces.
5. If you have a plated kit, select an index plate well from the Index well drop-down list and skip to step 8.
6. If applicable, specify an Index 1 sequence.

Select **Show Index Sequence/Show Index Names** to toggle between showing the name of the index and the index sequence.

7. If applicable, specify an Index 2 sequence.

i | During analysis, the iSeq™ 100, MiniSeq™, and NextSeq™ Systems and index kits automatically reverse complement the i5 indexes in custom library prep and index kits. Make sure that the i5 indexes are in the forward orientation.

8. [Optional] Enter a sample project name in the Sample Project field.
Use alphanumeric characters, dashes, underscores, or spaces.
9. [Optional] Select **Export Sample Sheet** to export sample information in *.csv format.
10. Select **Save Run**.

Import Sample Sheet

1. If you do not have a sample sheet to import, see [Enter Samples Manually on page 3](#) for instructions on how to create and export a sample sheet. Edit the file as follows.
 - a. Open the sample sheet in a text editor.
 - b. Enter the sample information in the [Data] section of the file.
 - c. Save the file. Make sure that the sample IDs are unique.
2. Select **Import Sample Sheet** at the top of the Create Run screen and browse to the location of the sample sheet.
Make sure that the information in the sample sheet is correct. Incorrect information can impact the sequencing run.
3. When finished, select **Save Run**.

Analysis Methods

The 16S Metagenomics analysis module performs the following analysis steps and then writes analysis output files to the folder.

- Demultiplexes index reads
- Generates FASTQ files

Demultiplexing

Demultiplexing compares each Index Read sequence to the index sequences specified for the run. No quality values are considered in this step.

Index reads are identified using the following steps:

- Samples are numbered starting from 1 based on the order they are listed for the run.
- Sample number 0 is reserved for clusters that were not assigned to a sample.

FASTQ File Generation

After demultiplexing, the software generates intermediate analysis files in the FASTQ format, which is a text format used to represent sequences. FASTQ files contain reads for each sample and the associated quality scores. Any controls used for the run and clusters that did not pass filter are excluded.

Each FASTQ file contains reads for only one sample, and the name of that sample is included in the FASTQ file name. FASTQ files are the primary input for alignment.

Classification of Reads

The classification step uses `ClassifyReads`, a proprietary algorithm that provides species-level classification for paired-end reads. This process involves matching short subsequences of the reads (called words) to a set of 16S reference sequences. The accumulated word matches for each read are used to assign reads to a particular taxonomic classification. Analysis results list the total number of classified clusters for each sample at each taxonomic level. Statistics are written to the file `Classification.txt`.

Current Taxonomy

The current taxonomy is stored in `Taxonomy.dat`. As of Local Run Manager v2.0 and v3.0, the 16S Metagenomics workflow generates classifications to the species level. For information about setting up analysis to genus level only, see [Sample Sheet Settings for Analysis on page 1](#).

The taxonomy database for the 16S Metagenomics workflow is an Illumina-curated version of the Greengenes database (greengenes.secondgenome.com/downloads/database/13_5). To generate species-level classifications, the following filters are applied:

- Filter all entries where the 16S sequence length was below 1250 bp.
- Filter all entries with more than 50 wobble bases (M, R, W, S, Y, K, V, H, D, B, and N).
- Filter all entries that are partially classified with no classification for genus or species.

The following taxonomic counts are available for the 16S Metagenomics workflow.

Taxonomy	Count
Kingdoms	3
Phyla	33
Classes	74
Orders	148
Families	321
Genera	1086
Species	6466

Alternative Taxonomy Database

You can prepare an alternative taxonomy database using the tool CreateTaxonomyDatabase distributed with Local Run Manager. This tool is in the Local Run Manager install folder, typically on the C: drive:

C:\Illumina\Local Run

Manager\Modules\MetagenomicsWorkflow\<>version>\Workflows\MetagenomicsWorker\CreateTaxonomyDatabase.exe.

CreateTaxonomyDatabase is a command-line tool; run it without arguments for a description of available options. For an example of a valid FASTA file, see:

greengenes.lbl.gov/Download/Sequence_Data/Fasta_data_files/current_GREENGENES_gg16S_unaligned.fasta.gz

The 16S Metagenomics workflow provides species-level classification. To configure the workflow for genus-level classification, use the TaxonomyFile sample sheet setting and specify `gg_13_5_genus_32bp.dat`. For more information, see [Sample Sheet Settings for Analysis on page 1](#).

Statistics Reporting

Statistics are summarized and reported, and written to the Alignment folder.

View Analysis Results

1. From the Local Run Manager dashboard, select the run name.
2. From the Run Overview tab, review the sequencing run metrics.
3. To change the analysis data file location for future requeues of the selected run, select the **Edit** icon, and edit the output run folder file path.
The file path leading up to the output run folder is editable. The output run folder name cannot be changed.
4. [Optional] Select the **Copy to Clipboard** icon to copy the output run folder file path.
5. Select the Sequencing Information tab to review run parameters and consumables information.
6. Select the Samples & Results tab to view the analysis report.
 - If analysis was requeued, select the appropriate analysis from the Select Analysis drop-down list.
7. [Optional] Select the **Copy to Clipboard** icon to copy the Analysis Folder file path.

Analysis Report

Analysis results are summarized on the Samples & Results tab. The report is also available in a PDF file format for each sample in the Analysis folder.

Sample Configuration

Row	Description
Sample ID	The sample ID from the sample sheet. Sample ID must always be a unique value.
Sample Name	The sample name from the sample sheet.
Run Folder	The folder location of the run.
Taxonomy File	Name of the taxonomy file in *.dat format.

Sample Information

Column	Description
Total Reads	The number of reads in a run.
Reads Passing Quality Filtering	The number of reads passing filter.
% Reads Passing Quality Filtering	The percentage of reads passing filter.

Classification Results

Report	Description
Classification Statistics	A table with the number of reads and percentage of reads classified to each classification level: Kingdom, Phylum, Class, Order, Family, Genus, Species.
Classification Rate by Taxonomic Level	A bar graph showing the percentage of reads classified to each classification level.

Report	Description
Sunburst Classification Chart	A sunburst chart showing the relative abundance of the classification results within each taxonomic level. Select on any category to zoom in, select the center button to zoom out, and select Reset Plot to return to a Kingdom-level view.
Top 20 Classification Results by Taxonomic Level	A bar chart showing the relative abundance of the top 20 classification results within each taxonomic level. Mouse over any category to see its description and abundance.
Classification Results by Taxonomic Level	Tables for each taxonomic level showing the number and percentage of reads for the highest eight taxonomic classifications. Pie charts for each taxonomic level show all classifications above 3.5% abundance. The "Other" category in the pie charts is the sum of all classifications with less than 3.50 % abundance.

Analysis Output Files

The following analysis output files are generated for the 16S Metagenomics analysis module and provide analysis results for alignment. Analysis output files are in the Alignment folder.

File Name	Description
Classification.txt	Contains the total number of classified clusters for each sample at each taxonomic level.
Demultiplexing (*.demux)	Intermediate files containing demultiplexing results.
DemuxSummary*	Intermediate files containing demultiplexing summary results.
FASTQ (*.fastq.gz)	Intermediate files containing quality-scored base calls. FASTQ files are the primary input for the alignment step.
*.rarefaction.txt	Contains the number of unique genera discovered by the number of reads classified.

Supplementary Output Files

The following output files provide supplementary information, or summarize run results and analysis errors. Although these files are not required for assessing analysis results, they can be used for troubleshooting purposes.

File Name	Description
AnalysisLog.txt	Processing log that describes every step that occurred during analysis of the current run folder. This file does not contain error messages.
AnalysisError.txt	Processing log that lists any errors that occurred during analysis. This file will be empty if no errors occurred.
Summary.txt	Contains summary statistics specific to the run.

Custom Analysis Settings

Custom analysis settings are intended for technically advanced users. If settings are applied incorrectly, serious problems can occur.

Add a Custom Analysis Setting

1. From the Module-Specific Settings section of the Create Run screen, select .
2. Select **+ Add custom setting**.
3. In the custom setting field, enter the setting name as listed in the Available Analysis Settings section.
4. In the setting value field, enter the setting value.
5. To remove a setting, select the delete icon .

Sample Sheet Settings for Analysis

Parameter	Description
Adapter	<p>Only specify for custom library prep and index kits.</p> <p>Specify the 5' portion of the adapter sequence to prevent reporting sequence beyond the sample DNA.</p> <p>Illumina recommends adapter trimming for Illumina DNA Prep (formally known as Nextera DNA), Nextera XT, and Nextera Mate Pair libraries.</p> <p>To specify two or more adapter sequences, separate the sequences by a plus (+) sign. For example:</p> <p>CTGTCTCTTATACACATCT+AGATGTGTATAAGAGACAG</p>
AdapterRead2	<p>Only specify for custom library prep kits and index kits.</p> <p>Specify the 5' portion of the Read 2-adapter sequence to prevent reporting sequence beyond the sample DNA.</p> <p>Use this setting to specify a different adapter other than the one specified in the Adapter setting.</p>
TaxonomyFile	<p>This setting overrides the taxonomy database; default is taxonomy.dat.</p> <p>Species-level classification is enabled, by default.</p>

Technical Assistance

For technical assistance, contact Illumina Technical Support.

Website: www.illumina.com
Email: techsupport@illumina.com

Illumina Technical Support Telephone Numbers

Region	Toll Free	International
Australia	+61 1800 775 688	
Austria	+43 800 006249	+43 1 9286540
Belgium	+32 800 77 160	+32 3 400 29 73
Canada	+1 800 809 4566	
China		+86 400 066 5835
Denmark	+45 80 82 01 83	+45 89 87 11 56
Finland	+358 800 918 363	+358 9 7479 0110
France	+33 8 05 10 21 93	+33 1 70 77 04 46
Germany	+49 800 101 4940	+49 89 3803 5677
Hong Kong, China	+852 800 960 230	
India	+91 8006500375	
Indonesia		0078036510048
Ireland	+353 1800 936608	+353 1 695 0506
Italy	+39 800 985513	+39 236003759
Japan	+81 0800 111 5011	
Malaysia	+60 1800 80 6789	
Netherlands	+31 800 022 2493	+31 20 713 2960
New Zealand	+64 800 451 650	
Norway	+47 800 16 836	+47 21 93 96 93
Philippines	+63 180016510798	
Singapore	1 800 5792 745	
South Korea	+82 80 234 5300	

Region	Toll Free	International
Spain	+34 800 300 143	+34 911 899 417
Sweden	+46 2 00883979	+46 8 50619671
Switzerland	+41 800 200 442	+41 56 580 00 00
Taiwan, China	+886 8 06651752	
Thailand	+66 1800 011 304	
United Kingdom	+44 800 012 6019	+44 20 7305 7197
United States	+1 800 809 4566	+1 858 202 4566
Vietnam	+84 1206 5263	

Safety data sheets (SDSs)—Available on the Illumina website at support.illumina.com/sds.html.

Product documentation—Available for download from support.illumina.com.



Illumina

5200 Illumina Way

San Diego, California 92122 U.S.A.

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America)

techsupport@illumina.com

www.illumina.com

For Research Use Only. Not for use in diagnostic procedures.

© 2022 Illumina, Inc. All rights reserved.

illumina[®]