

DRAGEN v3.3.5

Software Release Notes

April 19, 2019

Introduction

These release notes detail the key changes to software components for the Illumina® DRAGEN™ Bio-IT Platform since the package containing DRAGEN v3.2.8

If you are upgrading from a version prior to DRAGEN v3.2.8, please review the release notes for DRAGEN v3.2.8 for a list of features and bug fixes introduced in that version.

The software package includes:

- DRAGEN SW Intel Centos 6 - dragen-3.3.5.el6.x86_64
- DRAGEN SW Intel Centos 7 - dragen-3.3.5.el7.x86_64
- DRAGEN SW IBM Centos 7 - dragen-3.3.5.el7.ppc64le

Contents

Highlights	3
Fixes and Improvements	3
Summary of Changes	3
Somatic T/N Small Variant Calling	3
CNV Caller	3
Structural Variant Caller	4
Metrics	4
Repeat Expansion Calling	4
RNA Quantification	5
Known Limitations	5
SW Installation	5

Highlights

- Speed and Accuracy improvements to the Somatic T/N pipeline
 - Up to 6x speed improvement. Most 100x/40x T/N datasets now complete in under 2 hours.
 - Accuracy improvement compared to DRAGEN v3.2, on-par or better than GATK 4.1.
 - Accuracy has been validated against a range of tumor purities, library preps, and sequencing instruments.
- CNV DeNovo calling
 - Support for DeNovo calling and scoring with pedigree input.
- Speed improvement of BCL conversion
 - Up to 2x speed improvement on NovaSeq data when using DRAGEN Phase2 servers.
 - Up to 1.2-1.5x speed improvement on other servers.
- Structural Variants
 - DeNovo pedigree scoring.
 - Caller upgraded to Manta v1.5.1.
- Metrics
 - Added a model for detection of sample cross-contamination in human species.
 - Added MAPQ and BQ coverage filters for each selected coverage region.
- Repeat Expansion Calling
 - Updated repeat expansion caller to GraphExpansionHunter, which allows calling more complex repeat loci.
- Added RNA Quantification module to estimate transcript-level gene-expression results.

Fixes and Improvements

- Added support for quad and multi-generation pedigree calling in a single execution of the small variant joint caller.
- Added Beta support for read collapsing on the Illumina TSO-500 UMI design.
- Changes to the output VCFs and filenames are described below.

Summary of Changes

A summary of key changes is listed below. Please refer to the DRAGEN Bio-IT Platform User Guide for more information.

Somatic T/N Small Variant Calling

- Up to 6-fold speed improvement on datasets that were previously HMM-limited, with typical 100X/40X tumor-normal runs now finishing within 1h40m on a local server or 2h30m on AWS.
- Accuracy optimized for a broader range of datasets, with improvement for both snvs and indels on most datasets.

CNV Caller

- The CNV caller now supports de novo calling.
 - Multisample VCF support, starting from normalized signal files (*.tn.tsv) of single sample runs.
 - New *.tn.tsv files must be generated with this version of DRAGEN to be compatible with the de novo CNV caller.
 - De novo calling and scoring for valid trios defined in a pedigree file.
 - Multiple trios supported.

- Output VCF changes
 - The ID field in the output VCF now also encodes the contig of the event.
 - Now formatted as DRAGEN:<event>:<chr>:<start>-<stop>.
 - This is to comply with the VCF spec and ensure that the ID field is unique within the VCF.
 - Example: DRAGEN:LOSS:chr1:2841405-2847435
 - Added support for tabix of CNV VCFs.
- Filtering changes
 - Introduced a new option `cnv-filter-bin-support-ratio` to allow control of filtering events based on number of supporting bins.

Structural Variant Caller

- Structural Variant caller is updated to Manta 1.5.x
 - Improved accuracy, particularly improved precision for germline calling.
 - Improved runtime.
- Supports de novo calling
 - De novo scoring for a valid trio defined in a pedigree file.
 - Adds DQ and DN tags in the FORMAT field in the multi-sample VCF of germline calls.
- Output changes
 - VCF format changes
 - Change filters for easy interpretation of multi-sample germline variant vcf.
 - Add record-level filter 'SampleFT' when no sample passes all sample level filters.
 - Add sample-level filter 'HomRef' for homogyzous reference calls.
 - No more sample-level filters are applied at the record level even if it applies to all samples.
 - Change representation of inversions in the VCF output
 - Intrachromosomal translocations with inverted breakpoints are now reported as two breakend (BND) records.
 - Previously they were reported in the VCF using the inversion (INV) allele type.
 - The SV final output VCF is now available in the <output-directory>/<output-file-prefix>.sv.vcf.gz
 - The SV intermediate outputs moved to the <output-directory>/sv folder

Metrics

- The coverage report output file names have changed.
- All reports for `qc-coverage-region-i` are output in `qc-coverage-region-i*.bed` and `qc-coverage-region-i*.csv` files, where *i* can be 1, 2, or 3.

Repeat Expansion Calling

- Updated GraphExpansionHunter to allow calling more complex repeat loci, including loci with multiple flanking STRs and SNVs.
- The new version of repeat calling does not rely on unaligned (eg, fully in-repeat) reads. This makes calls more robust to sequencing bias, which means that repeat lengths will not be estimated beyond the library fragment size.
- This version uses a new repeat-spec (variant catalog) format; see the DRAGEN Bio-IT Platform User Guide and `repeat-specs/hg19/variant_catalog.json` for an example.
- Add SMA calling using the same graph alignment approach. This feature allows detecting absence of the fully functional allele at the duplicated SMN1/2 locus. Variant catalogs with SMN are found in `repeat-specs/experimental`.
- Realignment of repeat reads are now output in BAM format.

RNA Quantification

- Added support for quantification of gene expression from RNA-seq data.
- The module outputs the estimated expression of annotated transcripts and genes in Transcripts per Million (TPM) and read-counts units, using an EM method for deconvolution.
- Optionally includes GC-bias correction.
- Quantification can be enabled with map/align in RNA mode, by setting `--enable-rna-quantification` to true and supplying the transcript annotation file (GTF/GFF) with `--annotation-file`.
- Supports both stranded and un-stranded paired-end RNA-Seq protocol.
- RNA Quantification is still in Beta.

Known Limitations

- CNV de novo calling does not yet support Father/Proband calling on chrY.
- Structural Variant caller should be run with a BED file containing the set of regions to call, to avoid SV calls on alt and decoy contigs commonly found in hg38 references.

SW Installation

1. Install the appropriate release based on your Linux OS with the following command:

```
sudo sh <DRAGEN .run file>
```

2. Cold boot the server so that the new SW is fully installed with the updated FPGA HW image.

md5checksum:

```
70deab90191d781f3ad99a0758c1029b dragen-3.3.5.e16.x86_64.run  
d91b0a57ba24987006abd2f0db3376cb dragen-3.3.5.e17.x86_64.run  
ac249472668a8dbaf3dd1ccf59b4619d dragen-3.3.5.e17.ppc64le.run
```