

GenomeStudio v2.0.4

Release Notes

FOR RESEARCH USE ONLY

© 2017 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CPro, CytoChip, DesignStudio, Epicentre, GAIIX, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, Nextera, NextBio, NextSeq, NovaSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

Introduction

These Release Notes detail the key changes to the GenomeStudio Software since the release of GenomeStudio Software v2.0.3. Note that GenomeStudio Software v2.0.4 only includes the Genotyping Module and the Polyploid Genotyping Module. Users interested in running the Methylation Module and the Gene Expression Module should continue to run GenomeStudio Software v2011.1. GenomeStudio Software v2.0.3 and GenomeStudio Software v2011.1 can be installed and run simultaneously on the same system.

NEW FEATURES:

- None.

DEFECT REPAIRS:

- Removed the user authentication functionality. A recent update of the authentication platform is causing the authentication to fail preventing our user to complete the installation of the GenomeStudio software successfully.

KNOWN ISSUES:

- None