

Local Run Manager

Resequencing Analysis Module

Release Notes

3/14/2016

Introduction

These Release Notes detail notable items for the Local Run Manager Resequencing Analysis Module v1.1.0.9 release.

For more information about this analysis module and how to use it, refer to the Local Run Manager Resequencing Analysis Module Workflow Guide, available from the Local Run Manager Support Page, Documentation and Literature, on illumina.com.

http://support.illumina.com/sequencing/sequencing_software/local-run-manager/documentation.html

NEW FEATURES:

- None

DEFECT REPAIRS:

- Fixed issue causing all analyses using the Starling variant caller to abort

KNOWN ISSUES:

- None