

Local Run Manager

DNA Enrichment Analysis Module

Release Notes

11/FEB/2016

Introduction

These Release Notes detail notable items for the Local Run Manager DNA Enrichment Analysis Module v1.0 release.

For more information about this analysis module and how to use it, refer to the Local Run Manager DNA Enrichment Module Workflow Guide, available from the Local Run Manager Support Page, Documentation and Literature, on illumina.com.

http://support.illumina.com/sequencing/sequencing_software/local-run-manager/documentation.html

NEW FEATURES:

- Initial Release

DEFECT REPAIRS:

- None

KNOWN ISSUES:

- The Enrichment Sequencing Report html file incorrectly lists the software version numbers of Picard and GATK as "NA". The correct version numbers are as follows:

Software	Corrected Version Number
Picard (Calculate HS Metrics)	1.79(1282)
GATK (Variant Caller)	V1.6-23-gf0210b3