

Sequencing Analysis Viewer

Software Release Notes

SAV 2.5.12

November 2022

INTRODUCTION

These Release Notes detail the key changes to the Sequencing Analysis Viewer (SAV) Software since the release of SAV 2.4.7

Software Application	Prior Version	New Version
Sequencing Analysis Viewer	2.4.7	2.5.12

Sequencing Analysis Viewer version 2.5.12 is used to visualize sequencing run data from NextSeq 1000/2000. For data generated on other platforms please use the prior version.

This SAV version requires Windows. NET Framework 4.5.1 and Visual C++ redistribute for Visual Studio 2015. If your system is using Windows .NET Framework 4 or lower, continue using SAV 1.8.37.

NEW FEATURES

- Updated SAV to support NextSeq 1k2k long-read data
- This version of SAV is intended for NextSeq 1k2k customers with software v1.5 or greater. If you are using any other instruments, please use SAV v2.4.7